

# TYPHI<sub>LIZA</sub>



ELISA kit for the detection of specific IgM antibodies against a specific antigen for *Salmonella typhi*.

## FEATURES

- Convenient - can use EIA microtitre strip or well plate
- Ideal for private or hospital laboratories that screen large number of tests per day

## BENEFITS

- Early and definitive diagnosis of typhoid fever
- Results are easy to interpret
- Very small volume (1  $\mu$ l) of serum used
- High sensitivity and specificity
- Can screen a large number of samples at one time

## PERFORMANCE CHARACTERISTICS

Method	: Enzyme Immunoassay
Detection	: Specific human IgM antibodies against specific <i>Salmonella typhi</i> antigen in human serum
Specimen	: Serum
Specimen volume	: 1 $\mu$ l
Shelf life	: 12 months
Assay time	: 2.5 hours

## TEST PERFORMANCE

Sensitivity	: 97%
Specificity	: 96%
Positive Predictive Value	: 96%
Negative Predictive Value	: 97%

## APPLICATION

For the detection of acute typhoid fever in places such as hospitals or laboratories.

## ORDERING INFORMATION

MBDr/T04/96T



Reg. No. 06454 (Ukas)  
BS EN ISO 9001 : 2000 Certified

## PROCEDURE FOR TYPHI<sub>LIZA</sub>

- Prepare the wells in duplicate for Positive and Negative Controls and in single for Blank and Samples.
- Dilute samples in Diluent buffer (A1) to a final dilution of 1:100. Pipette 100  $\mu$ l of Controls and Samples into the corresponding wells. Pipette 100  $\mu$ l of Diluent buffer into the Blank well.
- Cover the tray with aluminium foil and incubate for 1 hour at 37°C.
- Wash the wells 4 times with 400  $\mu$ l of diluted Wash buffer. After the final washing cycle, turn the strips/plate down onto blotting paper or clean towel, and tap to remove any remaining solution.
- Add 100  $\mu$ l of Prediluted Anti-Human IgM\*HRP (B1) into each well.
- Cover the tray with aluminium foil and incubate for 1 hour at 37°C.
- Wash the wells as described in step (iv).
- Pipette 100  $\mu$ l of previously prepared ABTS Substrate Reagent into all wells. Cover the plate with aluminium foil and incubate at room temperature for 10 minutes, avoiding direct light exposure.
- After 10 minutes, add 100  $\mu$ l of Stop Solution (D1) into all well.
- Read the absorbance of the wells at 405 nm with a bichromatic spectrophotometer, preferably with reference wavelength at 492 nm (setting the instrument at zero with Blank well). Read the absorbance within 5 minutes after stopping the reaction.

## INTERPRETATION OF RESULTS

Test is valid if:

- Mean O.D of negative control is equal or lower than 0.200
- Mean O.D of positive control is equal or higher than 0.700

The cut-off value for the test is O.D. 0.300. An O.D. value above 0.300 is taken as positive ELISA result while an O.D. value below 0.300 is a negative result.