

# DOA<sup>SPOT TEST</sup>

## OPIATE<sup>SPOT TEST</sup>




MBDr's Opiate Spot Test is a rapid competitive binding immunoassay for qualitative determination of morphine, opiate and heroin in urine.

## THC<sup>SPOT TEST</sup>



MBDr's THC Spot Test is a rapid, qualitative, competitive binding immunoassay for the detection of tetrahydrocannabinol (THC) metabolites in urine.

## METHAMPHETAMINE<sup>SPOT TEST</sup>



MBDr's Methamphetamine Spot Test is a rapid qualitative, competitive binding immunoassay, for the determination of methamphetamine and its metabolites in urine.

## AMPHETAMINE<sup>SPOT TEST</sup>



MBDr's Amphetamine Spot Test is a rapid qualitative, competitive binding immunoassay, for the determination of amphetamine and its metabolites in urine.

### FEATURES

- No equipment needed
  - Built-in control
  - Easy to use and interpret
  - Results in 5 minutes
  - Room temperature storage
  - Ideal for field use
  - Detection sensitivity for opiate = 300 ng/ml or greater
  - Detection sensitivity for marijuana (THC) = 50 ng/ml or greater
  - Detection sensitivity for methamphetamine = 1000 ng/ml or greater
  - Detection sensitivity for amphetamine = 1000 ng/ml or greater
- [In accordance with the US Substance Abuse and Mental Health Services Administration (SAMHSA)]

### BENEFITS

- Simple and rapid
- No special equipment is needed
- Results are easy to interpret
- No sample preparation
- Results in 5 minutes

### PERFORMANCE CHARACTERISTICS

Method	: Immunochromatography
Detection	: Opiate : 300 ng/ml or >
	: Marijuana (THC) : 50 ng/ml or >
	: Methamphetamine : 1000 ng/ml or >
	: Amphetamine : 1000 ng/ml or >
Specimen	: Urine
Assay time	: 5 minutes
Reproducibility	: Very high
Precision	: Very good



### APPLICATION

Pre-employment or applicant testing  
 Scheduled testing  
 Random testing  
 Post-accident or for-cause testing  
 Treatment follow-up testing

### ORDERING INFORMATION

Opiate Spot Test	: MBD <sub>r</sub> /DOA01/50T
THC Spot Test	: MBD <sub>r</sub> /DOA02/50T
Methamphetamine Spot Test	: MBD <sub>r</sub> /DOA03/50T
Amphetamine Spot Test	: MBD <sub>r</sub> /DOA04/50T

### PROCEDURE FOR DOA SPOT TESTS (OPIATE, THC, METHAMPHETAMINE, AMPHETAMINE)



It is recommended that all positive results be confirmed by using an appropriate analytical technique

